


**Product Name: Recombinant Human TRIM5 Protein**  
**Catalog #: PEH3727**



## Summary

<b>Name</b>	TRIM5_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-MBP & C-His
<b>Accession #</b>	Q9C035
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	71 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	GHSFCQACTL ANHKKSM LDK GESSCPVCRI SYQPENIRPN RLVANIVEKL REVKLSPEGQ KVDHCA RHGE KLLFCQEDG KVICWLCERS QEHRGHHTFL TEEVAREYQV KLQAALEMLR QKQQEAELE ADIREEKASW KTQIQYDKTN VLADFEQLRD ILDWEE SNEL QNLEKEEDI LKSLTNSETE MVQQTQSLRE LISDLEHRLQ GVMELLQGV DGVIKRTENV TLKKPETFPK

## image



## Background

**Alternative Names** Tripartite motif-containing protein 5 (EC 2.3.2.27) (RING finger protein 88) (RING-type E3 ubiquitin transferase TRIM5)

## Note

**Product Name: Recombinant Human TRIM5 Protein**  
**Catalog #: PEH3727**



---

For Research Use Only , Not for Diagnostic Use