


Product Name: Recombinant Human TRIM5 Protein
Catalog #: PEH3728



Summary

Name	TRIM5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q9C035
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GHSFCQACTL ANHKKSM LDK GESSCPVCRI SYQPENIRPN RHVANIVEKL REVKLSPEGQ KVDHCA RHGE KLLFCQEDG KVICWLCERS QEHRGHHTFL TEEVAREYQV KLQAALEMLR QKQQEAELE ADIREEKASW KTQIQYDKTN VLADFEQLRD ILDWEE SNEL QNLEKEEDI LKSLTNSETE MVQQTQSLRE LISDLEHRLQ GVMELLQGV DGVIKRTENV TLKKPETFPK

image



Background

Alternative Names Tripartite motif-containing protein 5 (EC 2.3.2.27) (RING finger protein 88) (RING-type E3 ubiquitin transferase TRIM5)

Note

Product Name: Recombinant Human TRIM5 Protein
Catalog #: PEH3728



For Research Use Only , Not for Diagnostic Use