

Product Name: Recombinant Human CLSPN Protein
Catalog #: PEH3735



Summary

Name	CLSPN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HAW4
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DDDEEELLKR SEKLKQRMRL RKYLEDEAEV SGSDVGSSEDE YDGEEIDEYE EDVIDEVLPS DEELQSQIKK IHMKTMLDDD KRQLRLYQER YLADGDLHSD GPGRMRKFRW KNIDDASQMD LFHRDSDDDQ TEEQLDESEA RWRKERIERE QWLRDMAQQG KITAEDEEEI GEDSQFMILA KKVTAKALQK NASRPMVIQE SKSLLRNPFE AIRPGSAQQV KTGSLNQP AVLQKLAALS DHNPSAPRNS RNFVFHTLSP VKAEAAKES

image



Background

Alternative Names Claspin (hClaspin)

Product Name: Recombinant Human CLSPN Protein
Catalog #: PEH3735



Note

For Research Use Only , Not for Diagnostic Use