


**Product Name: Recombinant Human GOPC Protein**  
**Catalog #: PEH3742**



## Summary

<b>Name</b>	GOPC_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9HD26
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	21 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	QGVGPIRKVL LLKEDHEGLG ISITGGKEHG VPILISEIHP GQPADRCGGL HVGDAILAVN GVNLRDTKHK EAVTILSQQR GEIEFEVYV APEVDSDDEN VEYEDESGHR YRLYLDELEG GGNPGASCKD TSGEIKVLQG FNKKAVTDTH ENGLGTASE TPLDDGASKL DDLHTLYHKK

## image



## Background

<b>Alternative Names</b>	Golgi-associated PDZ and coiled-coil motif-containing protein (CFTR-associated ligand) (Fused in glioblastoma) (PDZ protein interacting specifically with TC10) (PIST)
--------------------------	--

## Note

**Product Name: Recombinant Human GOPC Protein**  
**Catalog #: PEH3742**



---

For Research Use Only , Not for Diagnostic Use