


Product Name: Recombinant Human BABAM1 Protein
Catalog #: PEH3753



Summary

Name	BABA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NWW8
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEEEEEEEHS AEPRPRTSN PEGAEDRAVG AQASVGSRS EGEAASADD GSLNTSGAGP KSWQVPPAP EVQIRTPRVN CPEKVIICLD LSEMSLPKL ESFNGSKTNA LNVSQKMIEM FVRTKHKIDK SHEFALVVVN DDTAWLSGLT SDPRELCACL

image



Background

Alternative Names	BRISC and BRCA1-A complex member 1 (Mediator of RAP80 interactions and targeting subunit of 40 kDa) (New component of the BRCA1-A complex)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human BABAM1 Protein
Catalog #: PEH3753

