


Product Name: Recombinant Human EIF2AK3 Protein
Catalog #: PEH3756



Summary

Name	E2AK3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NZJ5
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QEEDILLQR TQKTVRAVGP RSGNEKWNFS VGHFELRYIP DMETRAGFIE STFKPNENTE ESKIISDVEE QEAAIMDIVI KVSVDWKVM AFSKKGGHLE WEYQFCTPIA SAWLLKDGKV IPISLFDDTS YTSNDDVLED EEDIVEAARG ATENSVYLG M YRGQLYLQSS VRISEKFPSS PKALESVTNE

image



Background

Alternative Names	Eukaryotic translation initiation factor 2-alpha kinase 3 (EC 2.7.11.1) (PRKR-like endoplasmic reticulum kinase) (Pancreatic eIF2-alpha kinase) (HsPEK)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human EIF2AK3 Protein
Catalog #: PEH3756

