


Product Name: Recombinant Human BIN2 Protein
Catalog #: PEH3760



Summary

Name	BIN2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UBW5
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KLGKAVETKD ERFEQSASNF YQQQAEGHKL YKDLKNFLSA VKVMHESSKR VSETLQEIYS SEWDGHEELK AIVWNNDLLW EDYEEKLADQ AVRTMEIYVA QFSEIKERIA KRGRKLVLDYD SARHHLEAVQ NAKKKDEAKT AKAEFEFNKA QTVFEDLNQE LLEELPILYN SRIGCYVTIF QNISNLRDVF YREMSKLNHN LYEVMSKLEK QHSNKVFVVK GLSSSSRRSL VISPPVRTAT VSSPLTSPTS PSTLSLKSES ESVSATEDLA PDAAQGEDNS EIKELLEEEE IEKEGSEASS SEEDEPLPAC

image



Background

Alternative Names Bridging integrator 2 (Breast cancer-associated protein 1)

Product Name: Recombinant Human BIN2 Protein
Catalog #: PEH3760



Note

For Research Use Only , Not for Diagnostic Use