


Product Name: Recombinant Human IL20RA Protein
Catalog #: PEH3761



Summary

Name	I20RA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9UHF4
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PCVSGGLPKP ANITFLSINM KNVLQWTPPE GLQGVKVTYT VQYFIYGQKK WLNKSECRNI NRTYCDLSAE TSDYEHQYYA KVKAIWGTKC SKWAESGRFY PFLETQIGPP EVALTTDEKS ISVVLTAPEK WKRNPEDLPV SMQQIYSNLK YNVSVLNTKS NRTWSQCVTN HTLVLTWLEP NTLYCVHVES FVPGPPRRAQ PSEKQCARTL KDQSSEFKAK

image



Background

Alternative Names	Interleukin-20 receptor subunit alpha (IL-20 receptor subunit alpha) (IL-20R-alpha) (IL-20RA) (Cytokine receptor class-II member 8) (Cytokine receptor family 2 member 8) (CRF2-8) (IL-20R1) (ZcytoR7)
--------------------------	--

Product Name: Recombinant Human IL20RA Protein
Catalog #: PEH3761



Note

For Research Use Only , Not for Diagnostic Use