


Product Name: Recombinant Human IKZF2 Protein
Catalog #: PEH3766



Summary

Name	IKZF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UKS7
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ISSAYSQVYH PNRIERPISR ETADSHENNM DGPISLIRPK SRPQEREASP SNSCLDSTDS ESSHDDHQSY QGHPALNPKR KQSPAYMKED VKALDTTKAP KGSLKDIYKV FNGEGEQIRA

image



Background

Alternative Names Zinc finger protein Helios (Ikaros family zinc finger protein 2)

Note

For Research Use Only , Not for Diagnostic Use