


Summary

Name	MTUS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9ULD2
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QTAPDALPPE KTLELTQYKT KCENQSGFIL QLKQLLACGN TKFEALTVVI QHLLSEREEA LKQHKLSQE LVNLRGELVT ASTTCEKLEK ARNELQTVYE AFVQQHQA EK TERENRLKEF YTREYEKLRD TYIEEAKEYK MQLQE QFDNL NAAHETSKLE IEASHSEKLE LLKKAYEASL SEIKKGHEIE KKSLEDLLSE KQESLEKQIN DLKSENDALN EKLKSEEQKR RAREKANLKN PQIMYLEQEL ESLKAVLEIK NEKLHQQDIK LMKMEKLV DN NTALVDKLR FQQENEELKA RMDKHMAISR QLSTEQAVLQ ESLEKESKVN KRLSMENEEL LWKLHNGDLC SPKRSPTSSA

image



Background

Alternative Names	Microtubule-associated tumor suppressor 1 (AT2 receptor-binding protein) (Angiotensin-II type 2 receptor-interacting protein) (Mitochondrial tumor
--------------------------	---

Product Name: Recombinant Human MTUS1 Protein
Catalog #: PEH3767



suppressor 1)

Note

For Research Use Only , Not for Diagnostic Use