


Product Name: Recombinant Human EDAR Protein
Catalog #: PEH3768



Summary

Name	EDAR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UNE0
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EKDEFELTA TPAKPTKSEN DASSENEQLL SRSVDSDEEP APDKQGSPEL CLLSLVHLAR EKSATSNKSA GIQSRRKKIL DVYANVCGVV

image



Background

Alternative Names	Tumor necrosis factor receptor superfamily member EDAR (Anhidrotic ectodysplasin receptor 1) (Downless homolog) (EDA-A1 receptor) (Ectodermal dysplasia receptor) (Ectodysplasin-A receptor)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use