


Product Name: Recombinant Human EDAR Protein
Catalog #: PEH3769



Summary

Name	EDAR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UNE0
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EAPDNVVMFS EKDEFELTA TPAKPTKSEN DASSENEQLL SRSVDSDEEP APDKQGSPKL CLLSLVHLAR EKSATSNKSA GIQSRRKKIL DVYANVCGVV EGLSPTELPF

image



Background

Alternative Names	Tumor necrosis factor receptor superfamily member EDAR (Anhidrotic ectodysplasin receptor 1) (Downless homolog) (EDA-A1 receptor) (Ectodermal dysplasia receptor) (Ectodysplasin-A receptor)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human EDAR Protein
Catalog #: PEH3769

