


Product Name: Recombinant Human RBM19 Protein
Catalog #: PEH3775



Summary

Name	RBM19_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y4C8
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AAERSKTVIL VKNLPAGTLA AQLQETFGHF GSLGRVLLPE GGITAIVEFL EPLARKAFR HLAYSKFHHV PLYLEWAPVG VFSSTAPQKK KLQDTPSEPM EKDPAEPETV PDGETPEDEN PTEEGADNSS AKMEEEEEEE EEEEESLPGC TLFIKLNLFDT TTEEKLKEVF SKVGTVKSCS ISKKKNKAGV LLSMGFGFVE YRKPEQAQKA LKQLQGHVVD GHKLEVRIS

image



Background

Alternative Names	Probable RNA-binding protein 19 (RNA-binding motif protein 19)
--------------------------	--

Note

Product Name: Recombinant Human RBM19 Protein
Catalog #: PEH3775



For Research Use Only , Not for Diagnostic Use