


Product Name: Recombinant Human CD2AP Protein
Catalog #: PEH3779



Summary

Name	CD2AP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y5K6
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DAVHDELTI RVGEIIRNVK KLQEEGWLEG ELNGRRGMFP DNFVKEIKRE TEFKDDSLPI KRERHGNVAS LVQRISTYGL PAGGIQPHPQ TKNIKKTKK RQCKVLFYI PQNEDELELK VGDIIDINEE VEEGWWSGTL NNLGLFPSN FVKELEVTDG GETHEAQDDS ETVLAGPTSP IPSLGNVSET ASGSVTQPKK IRGIGFGDIF KEGSVKLRTR TSSSETEKK PEKPLILQSL GPKTQSVET KTDTEGKKA KEYCRTLFAY EGTNEDELTF KEGEIIHLIS KETGEAGWWR GELNGKEGVF PDNFAVQINE

image



Background

Alternative Names	CD2-associated protein (Adapter protein CMS) (Cas ligand with multiple SH3 domains)
--------------------------	---

Product Name: Recombinant Human CD2AP Protein
Catalog #: PEH3779



Note

For Research Use Only , Not for Diagnostic Use