


Product Name: Recombinant Human COLEC10 Protein
Catalog #: PEH3782



Summary

Name	COL10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y6Z7
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DSRPTAEVCA THTISPGPKG DDGEKGDPE EGKHGKVGGRM GPKGIKGELG DMGDQGNIGK TGPIGKKGDK GEKLLGIPG EKGKAGTVCD CGRYRKFVQ LDSIARLKT SMKFKNVIA GIRETEEFY YIVQEEKNYR ESLTHCRIRG GMLAMPKDEA ANTLIADYVA KSGFFRVFIG VNDLEREGQY MFTDNTPLQN YSNWNEGEPS DPYGHEDCVE MLSSGRWNDT ECHLTMVFVC

image



Background

Alternative Names Collectin-10 (Collectin liver protein 1) (CL-L1) (Collectin-34) (CL-34)

Note

Product Name: Recombinant Human COLEC10 Protein
Catalog #: PEH3782



For Research Use Only , Not for Diagnostic Use