


Product Name: Recombinant Human NFIC Protein
Catalog #: PEH3787



Summary

Name	NFIC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P08651
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RHKSGSMEED VDTSPGGDYY TSPSSPTSSS RNWTEDMEGG ISSPVKKTEM DKSPFNPSPP QDSPRLSSFT QHHRPVIAPH

image



Background

Alternative Names	Nuclear factor 1 C-type (NF1-C) (Nuclear factor 1/C) (CCAAT-box-binding transcription factor) (CTF) (Nuclear factor I/C) (NF-I/C) (NFI-C) (TGGCA-binding protein)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use