


Product Name: Recombinant Human NRCAM Protein
Catalog #: PEH3793



Summary

Name	NRCAM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92823
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KEKEDAHADP EIQPMKEDDG TFGEYSDAED HKPLKKGSRTPSDRTVKKED SDDSLVDYGE GVNGQFNEDG SFIGQYSGKK

image



Background

Alternative Names	Neuronal cell adhesion molecule (Nr-CAM) (Neuronal surface protein Bravo) (hBravo) (NgCAM-related cell adhesion molecule) (Ng-CAM-related)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use