



Product Name: Recombinant Human NECTIN1 Protein
Catalog #: PEH3816



Summary

Name	NECT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15223
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGTQAVLRAK KGQDDKVLVA TCTSANGKPP SVVSWETRLK GEAEYQEIRN PNGTVTVISR YRLVPSREAH QQSLACIVNY HMDRFKESLT LNVQYEPEVT IEGFDGNWYL QRMDVKLTCK ADANPPATEY HWTTLNGLSLP KGVEAQNRTL FFKGPINYSL AGTYICEATN PIGTRSGQVE VNITEFPYTP SPEHGRRAG

image



Background

Alternative Names	Nectin-1 (Herpes virus entry mediator C) (Herpesvirus entry mediator C) (HveC) (Herpesvirus Ig-like receptor) (HIgR) (Nectin cell adhesion molecule 1) (Poliovirus receptor-related protein 1) (CD antigen CD111)
--------------------------	---

Note

Product Name: Recombinant Human NECTIN1 Protein
Catalog #: PEH3816



For Research Use Only , Not for Diagnostic Use