


Summary

Name	CD244_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9BZW8
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SISGVPLQLQ P NSIQTKVDS IAWKKLLPSQ NGFHHLKWE NGLSPSNTSN DRFSFIVKNL SLLIKAAQQQ DSGLYCLEVT SIGKVTAT FQVVFESLL PDKVEKRLQ GQGKILDRGR CQVALSCLVS RDGNVSYAWY RGSKLIQTAG NLTYLDEEVD INGTHTYTCN VSNPVSWEH TLNLTQDCQN

image



Background

Alternative Names	Natural killer cell receptor 2B4 (NK cell activation-inducing ligand) (NAIL) (NK cell type I receptor protein 2B4) (NKR2B4) (h2B4) (SLAM family member 4) (SLAMF4) (Signaling lymphocytic activation molecule 4) (CD antigen CD244)
--------------------------	---

Note

Product Name: Recombinant Human CD244 Protein
Catalog #: PEH3818



For Research Use Only , Not for Diagnostic Use