


Product Name: Recombinant Human MSMB Protein
Catalog #: PEH3823



Summary

Name	MSMB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P08118
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PGDSTRKCMD LKGNKHPINS EWQTDNCETC TCYETEISCC TLVSTPVGVD KDNCQRIFKK EDCKYIVVEK

image



Background

Alternative Names	Beta-microseminoprotein (Immunoglobulin-binding factor) (IGBF) (PN44) (Prostate secreted seminal plasma protein) (Prostate secretory protein of 94 amino acids) (PSP-94) (PSP94) (Seminal plasma beta-inhibin)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use