


Product Name: Recombinant Human ITGA4 Protein
Catalog #: PEH3837



Summary

Name	ITA4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P13612
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LEEKQINCEV TDNSGVVQLD CSIGYIYVDH LSRIDISFLL DVSSLSRAEE DLSITVHATC ENEEEMDNLK

image



Background

Alternative Names Integrin alpha-4 (CD49 antigen-like family member D) (Integrin alpha-IV) (VLA-4 subunit alpha) (CD antigen CD49d)

Note

For Research Use Only , Not for Diagnostic Use