



Product Name: Recombinant Human PIK3AP1 Protein
Catalog #: PEH3838



Summary

Name	BCAP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZUJ8
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FLDFFPDWAH WQELTCDDEP ETYVAAVKKA ISEDSGCDV TDTEPEDEKV VSYKQQNLP TVTSPGNLMV VQPDRIRCGA ETTVYVIVRC KLDDRVATEA EFSPEDSPSV RMEAKVENEY TISVKAPNLS SGNVSLKIYS GDLVVCETVI SYYTDMEEIG NLLSNAANPV EFMCQAFKIV PYNTETLDKL LTESLKNNIP ASGLHLFGIN QLEEDMMTN QRDEELPTLL HFAAKYGLKN

image



Background

Alternative Names	Phosphoinositide 3-kinase adapter protein 1 (B-cell adapter for phosphoinositide 3-kinase) (B-cell phosphoinositide 3-kinase adapter protein 1)
--------------------------	---

Product Name: Recombinant Human PIK3AP1 Protein
Catalog #: PEH3838



Note

For Research Use Only , Not for Diagnostic Use