


Product Name: Recombinant Human FOXJ3 Protein
Catalog #: PEH3841



Summary

Name	FOXJ3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q9UPW0
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SDSPRSSLNN SLSDQSLASV NLNSVGSVHS YTPVTSHPEVSQSLTPQQQ PQYNLPERDK QLLFSEYNFE DLSASFRSLY KSVFEQSLSQ QGLMNIPSES SQQSHTSCTY QHSPSSTVST HPHSNQSSLS NSHGSGLNNTT GSNSVAQVSL SHPQMHTQPS PHPPHRPHGL PQHPQRSPPH APHPQQHSQQL QSPHPQHPSP HQHIQHHPNH QHQLTHQAP PPPQQVSCNS GVSNDWYATL DMLKESCRIA

image



Background

Alternative Names Forkhead box protein J3

Note

Product Name: Recombinant Human FOXJ3 Protein
Catalog #: PEH3841



For Research Use Only , Not for Diagnostic Use