


Product Name: Recombinant Human Tau Protein
Catalog #: PEH3846



Summary

Name	TAU_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P10636
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MEDHAGTYGL GDRKDQGGYT MHQDQEGDTD AGLKESPLQT PTEDGSEEPG SETSDAKSTP TAEDVTAPLV DEGAPGKQAA AQPHTIEPEG TTAEAGIGD TPSLEDEAAG HVTQARMVSK SKDGTGSDDK KAKGADGKTK IATPRGAAPP GQKGQANATR IPAKTPPAPK TPPSSGEPPK SGDRSGYSSP GSPGTPGSRS RTPSLTPPT REPKKVAVVR TPPKSPSSAK SRLQTAPVPM PDLKNVYSKI GSTENLKHQP GGGKVQIINK KLDLSNVQSK CGSKDNIHV PGGGSVQIVY KPVDSLKVTS KCGSLGNIHH KPGGGQVEVK SEKLDFKDRV QSKIGSLDNI

image



Background

Alternative Names	Microtubule-associated protein tau (Neurofibrillary tangle protein) (Paired helical filament-tau) (PHF-tau)
--------------------------	---

Product Name: Recombinant Human Tau Protein
Catalog #: PEH3846



Note

For Research Use Only , Not for Diagnostic Use