


Product Name: Recombinant Human STXBP1 Protein
Catalog #: PEH3848



Summary

Name	STXB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P61764
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QEVTRSLKDF SSSKRMNTGE KTTMRDLSQM LKKMPQYQKE LSKYSTHLHL AEDCMKHYQG TVDKLCRVEQ DLAMGTDAEG

image



Background

Alternative Names	Syntaxin-binding protein 1 (MUNC18-1) (N-Sec1) (Protein unc-18 homolog 1) (Unc18-1) (Protein unc-18 homolog A) (Unc-18A) (p67)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use