


Product Name: Recombinant Human ADGRV1 Protein
Catalog #: PEH3851



Summary

Name	AGRV1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WXG9
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TFGERCAQME PNALPFRGIY GISNLTWAVE EEDFEEQTLT LIFLDGERER KVSQILDDD EPEGQEFFYV FLTNPQGGAQ

image



Background

Alternative Names	Adhesion G-protein coupled receptor V1 (ADGRV1) (EC 3.4.-.-) (G-protein coupled receptor 98) (Monogenic audiogenic seizure susceptibility protein 1 homolog) (Usher syndrome type-2C protein) (Very large G-protein coupled receptor 1) [Cleaved into: ADGRV1 subunit alpha; ADGRV1 subunit beta]
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ADGRV1 Protein
Catalog #: PEH3851

