

Product Name: Recombinant Human LILRB2 Protein
Catalog #: PEH3858



Summary

Name	LIRB2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N423
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	WTSTQRKADF QHPAGAVGPE PTDRGLQWRS SPAADAQEEN LYAAVKDTQP EDGVEMDTRA AASEAPQDVT YAQLHSLTLR

image



Background

Alternative Names	Leukocyte immunoglobulin-like receptor subfamily B member 2 (LIR-2) (Leukocyte immunoglobulin-like receptor 2) (CD85 antigen-like family member D) (Immunoglobulin-like transcript 4) (ILT-4) (Monocyte/macrophage immunoglobulin-like receptor 10) (MIR-10) (CD antigen CD85d)
--------------------------	---

Note

Product Name: Recombinant Human LILRB2 Protein
Catalog #: PEH3858



For Research Use Only , Not for Diagnostic Use