


Product Name: Recombinant Human NOTCH3 Protein
Catalog #: PEH3864



Summary

Name	NOTC3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UM47
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SGRTPLHTAV TADAQGVFQI LIRNRSTDLD ARMADGSTAL ILAARLAVEG MVEELIASHA DVNAVDELGK SALHWAAAVN NVEATLALLK NGANKDMQDS KEETPLFLAA REGSYEAAKL

image



Background

Alternative Names Neurogenic locus notch homolog protein 3 (Notch 3) [Cleaved into: Notch 3 extracellular truncation; Notch 3 intracellular domain]

Note

For Research Use Only , Not for Diagnostic Use