

Product Name: Recombinant Human DRAXIN Protein
Catalog #: PEH3867



Summary

Name	DRAXI_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NBI3
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKEKHRGKLS SDGNETSPA E GPCHHHQDC LPGTCDDLRE HLCTPHNRGL NNKCFDDCMC VEGLRCYAKF HRNRRVTRRK GRCVEPETAN

image



Background

Alternative Names	Draxin (Dorsal inhibitory axon guidance protein) (Dorsal repulsive axon guidance protein) (Neucrin)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use