


Product Name: Recombinant Human AP3B1 Protein
Catalog #: PEH3870



Summary

Name	AP3B1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O00203
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QKPAPLLESP FKDRDHFQLG TLSHTLNIKA TGYLELSNWP EVAPDPSVRN VEVIELAKEW TPAGKAKQEN SAKKFYSESE EEEDSSDSSS

image



Background

Alternative Names AP-3 complex subunit beta-1 (Adaptor protein complex AP-3 subunit beta-1) (Adaptor-related protein complex 3 subunit beta-1) (Beta-3A-adaptin) (Clathrin assembly protein complex 3 beta-1 large chain)

Note

For Research Use Only , Not for Diagnostic Use