


Product Name: Recombinant Human CTSL Protein
Catalog #: PEH3873



Summary

Name	CATL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P07711
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GPQGNEGCNG GLMDYAFQYV QDNGGLDSEE SYPYEATEES CKYNPKYSVA NDTGFVDIPK QEKALMKAVA TVGPISVAID AGHESFLFYK EGIYFEPDCS SEDMDHGVLV

image



Background

Alternative Names	Procathepsin L (EC 3.4.22.15) (Cathepsin L1) (Major excreted protein) (MEP) [Cleaved into: Cathepsin L; Cathepsin L heavy chain; Cathepsin L light chain]
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use