


Product Name: Recombinant Human IL31RA Protein
Catalog #: PEH3884



Summary

Name	IL31R_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q8NI17
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VYYRKNLTC TWSPGKESY TQYTVKRTYA FGEKHDNCTT NSSTSENRAS CSFFLPRITI PDNYTIEVEA ENG DGVIKSH MTYWRLLENIA KTEPPKIFRV KPVLGIRMI QIEWIKPELA PVSSDLKYTL RFRTVNSTSW MEVNFANKRK DKNQTYNLTG LQPFTEYVIA LRCVAKESKF WSDWSQEKMG

image



Background

Alternative Names	Interleukin-31 receptor subunit alpha (IL-31 receptor subunit alpha) (IL-31R subunit alpha) (IL-31R-alpha) (IL-31RA) (Cytokine receptor-like 3) (GLM-R) (hGLM-R) (Gp130-like monocyte receptor) (Gp130-like receptor) (ZcytoR17)
--------------------------	--

Note

Product Name: Recombinant Human IL31RA Protein
Catalog #: PEH3884



For Research Use Only , Not for Diagnostic Use