


Product Name: Recombinant Human MAGED1 Protein
Catalog #: PEH3896



Summary

Name	MAGD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y5V3
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KSEMAFKAQN ATTKVGPNAT YNFSQSLNAN DLANSRPKTP FKAWNDDTTKA PTADTQTQNV NQAKMATSQA DIETDPGISE

image



Background

Alternative Names	Melanoma-associated antigen D1 (MAGE tumor antigen CCF) (MAGE-D1 antigen) (Neurotrophin receptor-interacting MAGE homolog)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use