


Product Name: Recombinant Human SMPD3 Protein
Catalog #: PEH3909



Summary

Name	NSMA2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NY59
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NHNQQDGDGSG SLGSPSASRE SLVKGRAGPD TSASGEPGAN SKLLYKASVV KAAAARRRRH PDEAFDHEVS AFFPANLDFL

image



Background

Alternative Names Sphingomyelin phosphodiesterase 3 (EC 3.1.4.12) (Neutral sphingomyelinase 2) (nSMase-2) (nSMase2) (Neutral sphingomyelinase II)

Note

For Research Use Only , Not for Diagnostic Use