


Product Name: Recombinant Human IBSP Protein
Catalog #: PEH3918



Summary

Name	SIAL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P21815
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GKGTSTTTTS PNGGFPTTP PQVYRTTSP FGKTTTVEYE GEYEYTGANE YDNGYEIYES ENGEPRGDNY RAYEDEYSYF KGQGYDGYDG

image



Background

Alternative Names	Bone sialoprotein 2 (Bone sialoprotein II) (BSP II) (Cell-binding sialoprotein) (Integrin-binding sialoprotein)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use