

Product Name: Recombinant Human CD6 Protein
Catalog #: PEH3920



Summary

| | |
|------------------------------|---|
| Name | CD6_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | N-SUMO & C-His |
| Accession # | P30203 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 53 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | PAPPDQLNTS SAESLWEPG ERLPVRLTNG SSSCSGTVEV RLEASWEPAC GALWDSRAAE AVCRALGCGG AEAASQLAPP TPELPPPPAA GNTSVAANAT LAGAPALLCS GAEWRLCEVV EHACRSDGRR ARVTCAENRA LRLVDGGGAC AGRVEMLEHG EWGSVCDDTW DLEDAHVVCR QLGCGWAVQA LPGLHFTPGR GPIHRDQVNC SGAEAYLWDC PGLPGQHYCG HKEDAGAVCS EHQSWRLTGG ADRCEGQVEV HFRGVWNTVC DSEWYPSEAK VLCQSLGCGT AVERPKGLPH SLSGRMYYS NGEELTSLNC SWRFNNSNLC SQSLAARVLC SASRSLHNLS TPEVPASVQT VTIESSVTVK |

image



Background

Alternative Names T-cell differentiation antigen CD6 (T12) (TP120) (CD antigen CD6) [Cleaved]

Product Name: Recombinant Human CD6 Protein
Catalog #: PEH3920



into: Soluble CD6]

Note

For Research Use Only , Not for Diagnostic Use