


Product Name: Recombinant Human IFNA8 Protein
Catalog #: PEH3926



Summary

Name	IFNA8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P32881
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FPQEEFDDKQ FQKAQAISVL HEMIQQTFNL FSTKDSSAAL DETLLDEFYI ELDQQLNDLE SCVMQEVGVI

image



Background

Alternative Names	Interferon alpha-8 (IFN-alpha-8) (Interferon alpha-B) (LeIF B) (Interferon alpha-B2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use