


Product Name: Recombinant Human CD83 Protein
Catalog #: PEH3927



Summary

Name	CD83_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q01151
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PEVKVACSED VDLPCTAPWD PQVPYTVSWV KLEGGGERM ETPQEDHLRG QHYHQKQQNG SFDAPNERPY SLKIRNTTSC NSGTYRCTLQ

image



Background

Alternative Names CD83 antigen (hCD83) (B-cell activation protein) (Cell surface protein HB15) (CD antigen CD83)

Note

For Research Use Only , Not for Diagnostic Use