



## Summary

<b>Name</b>	ZBT7A_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	O95365
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	16 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	AEGEDGDGPD VDGLAASTLL QQMMSSVGRA GAAAGDSDEE SRADDKGVMMD YYLKYFSGAH DGDVYPAWSQ KVEKKIRAKA FQKCPICEKV IQGAGKLPKH IRHTHTGEKPY ECNICKVRFT RQDKLKVHMR

## image



## Background

<b>Alternative Names</b>	Zinc finger and BTB domain-containing protein 7A (Factor binding IST protein 1) (FBI-1) (Factor that binds to inducer of short transcripts protein 1) (HIV-1 1st-binding protein 1) (Leukemia/lymphoma-related factor) (POZ and Krueppel erythroid myeloid ontogenic factor) (POK erythroid myeloid ontogenic factor) (Pokemon) (Pokemon 1) (TTF-I-interacting peptide 21) (TIP21) (Zinc finger protein 857A)
--------------------------	---

**Product Name: Recombinant Human ZBTB7A Protein**  
**Catalog #: PEH3933**

---



**Note**

For Research Use Only , Not for Diagnostic Use