


Product Name: Recombinant Human TMPO Protein
Catalog #: PEH3935



Summary

| | |
|------------------------------|--|
| Name | LAP2B_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | P42167 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 12 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | GPPDFSSDEE REPTVLGSG AAAAGRSRAA VGRKATKKT D KPRQEDKDDL DVTELTNEDL LDQLVKYGVN PGPIVGTTRK LYEKLLKLR EQGTESRSST |

image



Background

| | |
|--------------------------|--|
| Alternative Names | Lamina-associated polypeptide 2, isoforms beta/gamma (Thymopietin, isoforms beta/gamma) (TP beta/gamma) (Thymopietin-related peptide isoforms beta/gamma) (TPRP isoforms beta/gamma) [Cleaved into: Thymopietin (TP) (Splenin); Thymopentin (TP5)] |
|--------------------------|--|

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TMPO Protein
Catalog #: PEH3935

