


Product Name: Recombinant Human SLC1A4 Protein
Catalog #: PEH3936



Summary

| | |
|------------------------------|--|
| Name | SATT_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | P43007 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 10 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | VDWIVDRTTT VVNVEGDALG AGILHHLNQK ATKKGEQELA EVKVEAIPNC KSEETSPLV THQNPAGPVA SAPELESKES |

image



Background

Alternative Names Neutral amino acid transporter A (Alanine/serine/cysteine/threonine transporter 1) (ASCT-1) (SATT) (Solute carrier family 1 member 4)

Note

For Research Use Only , Not for Diagnostic Use