


Product Name: Recombinant Human DCBLD2 Protein
Catalog #: PEH3939



Summary

Name	DCBD2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96PD2
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GQENSWKPKK ARLKKPGPPW AAFATDEYQW LQIDLNKEKK ITGIITGST MVEHNYVSA YRILYSDDGQ KWTVYREPGV EQDKIFQGNK DYHQDVRNNF LPPIIARFIR VNPTQWQQKI AMKMELLGCQ FIPKGRPPKL TQPPPPRNSN DLKNTTAPPK IAKGRAPKFT QPLQPRSSNE FPAQTEQTTA SPDIRNTTVT

image



Background

Alternative Names	Discoidin, CUB and LCCL domain-containing protein 2 (CUB, LCCL and coagulation factor V/VIII-homology domains protein 1) (Endothelial and smooth muscle cell-derived neuropilin-like protein)
--------------------------	---

Note

Product Name: Recombinant Human DCBLD2 Protein
Catalog #: PEH3939



For Research Use Only , Not for Diagnostic Use