


Product Name: Recombinant Human TDRKH Protein
Catalog #: PEH3941



Summary

Name	TDRKH_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q9Y2W6
Host	E.coli
Species	Human
Molecular Mass	71 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	REERLTFVGE DDIEIEMRVP QEAVKLIIGR QGANIKQLRK QTGARIDVDT EDVGDERVLL ISGFPVQVCK AKAAIHQILT ENTPVSEQLS VPQRSVGRRI GRGGETIRSI CKASGAKITC DKESEGTLLL SRLIKISGTQ KEVAAAKHLI LEKVSEDEEL RKRIAHSAET RVPRKQPISV RREDMTEPGG AGEPALWKNT SSSMEPTAPL VTPPPKGGGD MAVVVSKEGS WEKPSDDSFQ KSEAQAIPEM

image



Background

Alternative Names	Tudor and KH domain-containing protein (Tudor domain-containing protein 2)
--------------------------	--

Note

Product Name: Recombinant Human TDRKH Protein
Catalog #: PEH3941



For Research Use Only , Not for Diagnostic Use