


Product Name: Recombinant Human CDCP1 Protein
Catalog #: PEH3952



Summary

Name	CDCP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H5V8
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GPAVGIYNDN INTEMPRQPK KFKQGRKDND SHVYAVIEDT MVYGHLLQDS SGSFLQPEVD TYRPFQGTMG VCPPSPPTIC SRAPTAKLAT EEPSPRPPE SESEPYTFESH PNNGDVSSKD

image



Background

Alternative Names	CUB domain-containing protein 1 (Membrane glycoprotein gp140) (Subtractive immunization M plus HEp3-associated 135 kDa protein) (SIMA135) (Transmembrane and associated with src kinases) (CD antigen CD318)
--------------------------	--

Note

Product Name: Recombinant Human CDCP1 Protein
Catalog #: PEH3952



For Research Use Only , Not for Diagnostic Use