


Product Name: Recombinant Human ALDH3A1 Protein
Catalog #: PEH3962



Summary

Name	AL3A1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P30838
Host	E.coli
Species	Human
Molecular Mass	59 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RPLQFRIQQL EALQRLIQEQ EQELVGALAA DLHKNEWNAY YEEVVYVLEE IEYMIQKLPE WAADEPVEKT PQTQQDELYI HSEPLGVVLV IGTWNYPFNL TIQPMVGAIA AGNSVVLKPS ELSENMASLL ATIIPQYLDK

image



Background

Alternative Names	Aldehyde dehydrogenase, dimeric NADP-preferring (EC 1.2.1.5) (ALDHIII) (Aldehyde dehydrogenase 3) (Aldehyde dehydrogenase family 3 member A1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use