


Product Name: Recombinant Human LILRB1 Protein
Catalog #: PEH3970



Summary

Name	LIRB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NHL6
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HWTSTQRKAD FQHPAGAVGP EPTDRGLQWR SSPAADAQEE NLYAAVKHTQ PEDGVEMDTR SPHDEDPQAV TYAEVKHSRP RREMASPPSP LSGEFLDTKD

image



Background

Alternative Names	Leukocyte immunoglobulin-like receptor subfamily B member 1 (LIR-1) (Leukocyte immunoglobulin-like receptor 1) (CD85 antigen-like family member J) (Immunoglobulin-like transcript 2) (ILT-2) (Monocyte/macrophage immunoglobulin-like receptor 7) (MIR-7) (CD antigen CD85j)
--------------------------	---

Note

Product Name: Recombinant Human LILRB1 Protein
Catalog #: PEH3970



For Research Use Only , Not for Diagnostic Use