


Product Name: Recombinant Human RNF43 Protein
Catalog #: PEH3991



Summary

Name	RNF43_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q68DV7
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NCTDISLQGV HGSSSTFCSS LSSDFDPLVY CSPKGDQQRV DMQPSVTSRP RSLDSVVPTG ETQVSSHVHY HRHRHHYK RFQWHGRKPG PETGVPQSRP PIPRTQPQPE PPSPDQQVTR

image



Background

Alternative Names E3 ubiquitin-protein ligase RNF43 (EC 2.3.2.27) (RING finger protein 43) (RING-type E3 ubiquitin transferase RNF43)

Note

For Research Use Only , Not for Diagnostic Use