



Product Name: Recombinant Human C5 Protein
Catalog #: PEH3994



Summary

Name	CO5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P01031
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QLVGGVPVTL NAQTIDVNQE TSDLDPKSV TRVDDGVASF VLNLPSGVTV LEFNVKTDAP DLPEENQARE GYRAIAYSSL SQSYLYIDWT

image



Background

Alternative Names	Complement C5 (C3 and PZP-like alpha-2-macroglobulin domain-containing protein 4) [Cleaved into: Complement C5 beta chain; Complement C5 alpha chain; C5a anaphylatoxin; Complement C5 alpha' chain]
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use