


Product Name: Recombinant Human PI16 Protein
Catalog #: PEH3997



Summary

Name	PI16_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6UXB8
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KVTDKTKVPS RSPENSLDPK MSLTGARELL PHAQEEAEAE AELPPSSEVL ASVFPAQDKP GELQATLDHT GHTSSKSLPN

image



Background

Alternative Names	Peptidase inhibitor 16 (PI-16) (Cysteine-rich secretory protein 9) (CRISP-9) (PSP94-binding protein) (CD antigen CD364)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use