


Product Name: Recombinant Human CHST3 Protein
Catalog #: PEH3998



Summary

Name	CHST3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7LGC8
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KLKQIPQALA DANSTDPALI LAENASLLSL SELDSAFSQL QSRLRNLSLQ LGVEPAMEAA GEEEEEQRKE

image



Background

Alternative Names	Carbohydrate sulfotransferase 3 (EC 2.8.2.17) (Chondroitin 6-O-sulfotransferase 1) (C6ST-1) (Chondroitin 6-sulfotransferase) (Galactose/N-acetylglucosamine/N-acetylglucosamine 6-O-sulfotransferase 0) (GST-0)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use